

## Anti-Myogenin/MYOG Antibody Picoband® HRP Conjugated

Catalog Number: A03238-1-HRP

### About MYOG

Myogenin, is a transcriptional activator encoded by the MYOG gene. Myogenin is a muscle-specific transcription factor that can induce myogenesis in a variety of cell types in tissue culture. It is a member of a large family of proteins related by sequence homology, the helix-loop-helix (HLH) proteins. It is essential for the development of functional skeletal muscle.

### Overview

Product Name	Anti-Myogenin/MYOG Antibody Picoband® HRP Conjugated
Reactive Species	Monkey, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P15173

### Technical Details

Immunogen	E.coli-derived human Myogenin/MYOG recombinant protein (Position: M1-N224).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Myogenin/MYOG Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.