

## Anti-CaMKII Alpha / Delta (phospho-T286) CAMK2A Antibody

Catalog Number: A03241T286

### Overview


Product Name	Anti-CaMKII Alpha / Delta (phospho-T286) CAMK2A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CaMKII Alpha / Delta (phospho-T286) CAMK2A Antibody catalog # A03241T286. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UQM7/Q13557

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human CaMKIIalpha around the phosphorylation site of Threonine 286.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-CaMKII Alpha / Delta (phospho-T286) CAMK2A Antibody (A03241T286) Images

---



Immunohistochemistry (IHC) analyzes of p-CaMKIIalpha/delta (T286) pAb in paraffin-embedded human brain tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CaMKII Alpha / Delta (phospho-T286) CAMK2A Antibody

For Research Use Only. Not for use in diagnostic procedures.