

Anti-LEKTI (H525) SPINK5 Antibody

Catalog Number: A03254-1


Overview

Product Name	Anti-LEKTI (H525) SPINK5 Antibody
Reactive Species	Human
Description	Boster Bio Anti-LEKTI (H525) SPINK5 Antibody catalog # A03254-1. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NQ38

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 496-543 of Human LEKTI.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-LEKTI (H525) SPINK5 Antibody (A03254-1) Images



Immunohistochemistry (IHC) analyzes of LEKTI (H525) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LEKTI (H525) SPINK5 Antibody

For Research Use Only. Not for use in diagnostic procedures.