

Anti-Pericentrin PCNT Antibody

Catalog Number: A03256-1

Overview

Product Name	Anti-Pericentrin PCNT Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Pericentrin PCNT Antibody catalog # A03256-1. Tested in IHC applications. This antibody reacts with Human, Mouse.
Application	IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95613

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	PCNT antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC-p 1:50-300

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

