

Anti-CD94 KLRD1 Monoclonal Antibody

Catalog Number: A03257-1


Overview

Product Name	Anti-CD94 KLRD1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-CD94 KLRD1 Monoclonal Antibody catalog # A03257-1. Tested in WB, Flow Cytometry, ELISA applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, WB
Clonality	Monoclonal 1H1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q13241


Technical Details

Immunogen	Purified recombinant fragment of human CD94 expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000

Anti-CD94 KLRD1 Monoclonal Antibody (A03257-1) Images



Western Blot analysis using CD94 Monoclonal Antibody against recombinant protein.



Flow cytometric analysis of RAJI cells using CD94 Monoclonal Antibody (green) and negative control (purple).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD94 KLRD1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.