

## Anti-Claudin 5 CLDN5 Antibody

Catalog Number: A03260-1


### Overview

Product Name	Anti-Claudin 5 CLDN5 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Claudin 5 CLDN5 Antibody catalog # A03260-1. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00501


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:2,000-1:5,000 IHC: 1:50-1:200

## Anti-Claudin 5 CLDN5 Antibody (A03260-1) Images



Western blot analysis of Claudin 5 on human lung cells lysates using anti-Claudin 5 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-Claudin 5 antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Claudin 5 CLDN5 Antibody

For Research Use Only. Not for use in diagnostic procedures.