

## Anti-PHF8 Monoclonal Antibody

Catalog Number: A03288

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-PHF8 Monoclonal Antibody  |
| Reactive Species     | Human, Mouse, Pig, Rat   |
| Description          | Boster Bio Anti-PHF8 Monoclonal Antibody catalog # A03288. Tested in WB applications. This antibody reacts with Human, Mouse, Pig, Rat.    |
| Application          | WB   |
| Clonality            | Monoclonal   |
| Formulation          | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.               |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Mouse  |
| Uniprot ID           | Q9UPP1   |

### Technical Details

|                     |                    |
|---------------------|--------------------|
| Isotype             | IgG                |
| Form                | Liquid             |
| Concentration       | 1 mg/mL            |
| Purification        | Protein A purified |
| Suggested Dilutions | WB, 1:1000-1:2000  |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.