

Anti-PHF8 Monoclonal Antibody

Catalog Number: A03288


Overview

Product Name	Anti-PHF8 Monoclonal Antibody
Reactive Species	Human, Mouse, Pig, Rat
Description	Boster Bio Anti-PHF8 Monoclonal Antibody catalog # A03288. Tested in WB applications. This antibody reacts with Human, Mouse, Pig, Rat.
Application	WB
Clonality	Monoclonal
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9UPP1

Technical Details

Immunogen	Purified recombinant human PHF8 (C-terminal) protein fragments expressed in E.coli.
Predicted Reactive Species	Canine
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL.
Purification	Protein A purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB, 1:1000-1:2000</p>

Anti-PHF8 Monoclonal Antibody (A03288) Images



Western Blot (WB) analysis using PHF8 Monoclonal antibody against A451,293 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PHF8 Monoclonal Antibody