

Anti-Calpain 10 CAPN10 Antibody

Catalog Number: A03289


Overview

Product Name	Anti-Calpain 10 CAPN10 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Calpain 10 CAPN10 Antibody catalog # A03289. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9HC96

Technical Details

Immunogen	Synthesized peptide derived from the N-terminal region of human Calpain 10.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL.
Purification	Immunogen affinity purified
Suggested Dilutions	WB, 1:500-1:2000

Anti-Calpain 10 CAPN10 Antibody (A03289) Images



Western Blot (WB) analysis of LOVO cells using Calpain 10 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calpain 10 CAPN10 Antibody