

## Anti-MOG Antibody

Catalog Number: A03294-1

### Overview


Product Name	Anti-MOG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MOG Antibody catalog # A03294-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16653

### Technical Details


Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

## Anti-MOG Antibody (A03294-1) Images

---



Immunohistochemical analysis of paraffin-embedded rat spinal cord tissue using anti-Myelin oligodendrocyte glycoprotein antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-Myelin oligodendrocyte glycoprotein antibody. Counter stained with hematoxylin.

---

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MOG Antibody

For Research Use Only. Not for use in diagnostic procedures.