

Anti-RASGRF1 Antibody

Catalog Number: A03338-2


Overview

Product Name	Anti-RASGRF1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RASGRF1 Antibody catalog # A03338-2. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13972


Technical Details

Immunogen	Recombinant full length Human RASGRF1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-RASGRF1 Antibody (A03338-2) Images



Western blot analysis of extract of Mouse brain, using RASGRF1 antibody.



Immunofluorescence analysis of U20S cell using RASGRF1 antibody. Blue: DAPI for nuclear staining.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RASGRF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.