

## Anti-Steroid hormone receptor ERR1 ESRRRA Antibody

Catalog Number: A03339-1


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-Steroid hormone receptor ERR1 ESRRRA Antibody  |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio Anti-Steroid hormone receptor ERR1 ESRRRA Antibody catalog # A03339-1. Tested in WB,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | Flow Cytometry, IP, WB  |
| Clonality            | Polyclonal  |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                |
| Host                 | Rabbit  |
| Uniprot ID           | P11474  |


### Technical Details

|                     |                                       |
|---------------------|---------------------------------------|
| Immunogen           | recombinant protein                   |
| Isotype             | IgG                                   |
| Form                | Liquid                                |
| Concentration       | 1 mg/ml                               |
| Purification        | ProA affinity purified                |
| Suggested Dilutions | WB: 1:1,000-1:2,000<br>FC: 1:50-1:100 |

## Anti-Steroid hormone receptor ERR1 ESRRRA Antibody (A03339-1) Images



Western blot analysis of Estrogen Related Receptor alpha on MCF-7 cell lysates using anti-Estrogen Related Receptor alpha antibody at 1/1,000 dilution.



Flow cytometric analysis of Hela cells with Estrogen Related Receptor alpha antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Steroid hormone receptor ERR1 ESRRRA Antibody

For Research Use Only. Not for use in diagnostic procedures.