

Anti-GlyT2/SLC6A5 Antibody

Catalog Number: A03349-1

About SLC6A5

GlyT2 is a transporter protein of the SLC6 family of sodium- and chloride-dependent neurotransmitter transporters. Like other members of the family, GlyT2 is trafficked to the cell surface as homo-oligomers from the endoplasmic reticulum (ER) (PMID: 32714574). GlyT2 is expressed in medulla, and to a lesser extent in spinal cord and cerebellum. GlyT2 has 3 isoforms with the molecular mass of 63, 86 and 87 kDa. The 63 kDa band is the unglycosylated GlyT2 polypeptide, whereas the 70-110 kDa band is a differential glycosylation form of GlyT2 (PMID: 17851090, 17962467).

Overview

Product Name	Anti-GlyT2/SLC6A5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GlyT2/SLC6A5 Antibody catalog # A03349-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y345

Technical Details

Immunogen	E.coli-derived human GlyT2/SLC6A5 recombinant protein (Position: C3-K197).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GlyT2/SLC6A5 Antibody

For Research Use Only. Not for use in diagnostic procedures.