

Anti-Paired box protein Pax-9 PAX9 Antibody

Catalog Number: A03356


Overview

| | |
|----------------------|--|
| Product Name | Anti-Paired box protein Pax-9 PAX9 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Paired box protein Pax-9 PAX9 Antibody catalog # A03356. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P55771 |

Technical Details

| | |
|---------------------|------------------------|
| Immunogen | recombinant protein |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | ProA affinity purified |
| Suggested Dilutions | WB: 1:1,000-1:2,000 |

Anti-Paired box protein Pax-9 PAX9 Antibody (A03356) Images



Western blot analysis of PAX9 on different lysates using anti-PAX9 antibody at 1/1,000 dilution. Positive control: Lane 1: HepG2 Lane 2: HeLa

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Paired box protein Pax-9 PAX9 Antibody

For Research Use Only. Not for use in diagnostic procedures.