

Anti-Phospho-SHIP-1 (Y1021) INPP5D Antibody

Catalog Number: A03358Y1021

Overview

Product Name	Anti-Phospho-SHIP-1 (Y1021) INPP5D Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-SHIP-1 (Y1021) INPP5D Antibody catalog # A03358Y1021. Tested in WB, IF, ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92835

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human SHIP1 around the phosphorylation site of Tyr1021. AA range:987-1036
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 IF/ICC 1:50-200 ELISA 1:2000-20000 IHC 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-SHIP-1 (Y1021) INPP5D Antibody

For Research Use Only. Not for use in diagnostic procedures.