

Anti-RNA demethylase ALKBH5 ALKBH5 Antibody

Catalog Number: A03360-2

Overview

Product Name	Anti-RNA demethylase ALKBH5 ALKBH5 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RNA demethylase ALKBH5 ALKBH5 Antibody catalog # A03360-2. Tested in IHC, WB applications. This antibody reacts with Human, Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6P6C2

Technical Details

Immunogen	Synthesized peptide derived from human ALKBH5
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used IHC-p 1:50-200, WB 1:500-2000 For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RNA demethylase ALKBH5 ALKBH5 Antibody