

Anti-Autophagy-related protein 9B ATG9B Antibody

Catalog Number: A03381-1


Overview

Product Name	Anti-Autophagy-related protein 9B ATG9B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Autophagy-related protein 9B ATG9B Antibody catalog # A03381-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q674R7


Technical Details

Immunogen	Recombinant fusion protein of human ATG9B(NP_775952.4).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200


Anti-Autophagy-related protein 9B ATG9B Antibody (A03381-1) Images



Western blot analysis of extracts of HUH-7 cells, using ATG9B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.



Immunofluorescence analysis of human placenta cells using ATG9B Rabbit pAb at dilution of 1:50 . Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using ATG9B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 30s.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Autophagy-related protein 9B ATG9B Antibody

For Research Use Only. Not for use in diagnostic procedures.