

## Anti-TRPM6 (Chak2) Antibody

Catalog Number: A03390

### Overview

Product Name	Anti-TRPM6 (Chak2) Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-TRPM6 (Chak2) Antibody (Catalog# A03390). Tested in IHC applications. This antibody reacts with Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Storage Instructions	Store at 4°C for short term. Aliquot and store at -20°C for long term. Avoid repeated freeze-thaw cycles.
Host	Guinea Pig
Uniprot ID	Q9BX84

### Technical Details

Immunogen	Mouse TRPM6-C Terminus Peptide-NH <sub>2</sub> - (GC)ERDKNRSSLERHTRL-OH
Isotype	Guinea Pig IgG
Form	Liquid
Purification	Whole serum
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

## Anti-TRPM6 (Chak2) Antibody (A03390) Images




Figure 1. Immunofluorescent staining data of TRPM6 using Anti-TRPM6 (Chak2) Antibody (A03390).

TRPM6 staining cryostat section of rat kidney. Tissues were stained with using Cy3-conjugated secondary antibodies.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRPM6 (Chak2) Antibody