

Anti-NG2 CSPG4 Antibody

Catalog Number: A03394-1


Overview

Product Name	Anti-NG2 CSPG4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-NG2 CSPG4 Antibody catalog # A03394-1. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6UVK1


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 FC: 1:50-1:100 IHC: 1:50-1:200

Anti-NG2 CSPG4 Antibody (A03394-1) Images



Immunohistochemical analysis of paraffin-embedded rat brain tissue using anti-NG2 antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti-NG2 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NG2 CSPG4 Antibody

For Research Use Only. Not for use in diagnostic procedures.