

Anti-IL12RB1 Antibody

Catalog Number: A03401


Overview

Product Name	Anti-IL12RB1 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-IL12RB1 Antibody catalog # A03401. Tested in WB,IF,IHC,ELISA applications. This antibody reacts with Mouse,Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P42701

Technical Details

Immunogen	Recombinant full length protein of human CD212.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:100-1:500 IF: 1:100-1:500 ELISA: 1:5000-1:10000

Anti-IL12RB1 Antibody (A03401) Images



Western blot analysis of CD212 expression in A431 (A), MCF7 (B), mouse liver (C), mouse thymus (D) whole cell lysates.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL12RB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.