

Anti-TCTP TPT1 Antibody

Catalog Number: A03414-1


Overview

Product Name	Anti-TCTP TPT1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TCTP TPT1 Antibody catalog # A03414-1. Tested in WB,ICC/IF,IHC,Flow Cytometry,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P13693


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC/IF: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-TCTP TPT1 Antibody (A03414-1) Images



Western blot analysis of TCTP on different lysates using anti-TCTP antibody at 1/1,000 dilution. Positive control: Lane 1: Mouse placenta tissue Lane 2: Mouse testis tissue Lane 3: HeLa Lane 4: A431 Lane 5: MCF-7



ICC staining TCTP in A431 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TCTP TPT1 Antibody

For Research Use Only. Not for use in diagnostic procedures.