

Anti-TIEG-1/2 KLF10 Antibody

Catalog Number: A03419-1


Overview

Product Name	Anti-TIEG-1/2 KLF10 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TIEG-1/2 KLF10 Antibody catalog # A03419-1. Tested in ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13118

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human KLF10/11. AA range:391-440
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 ELISA 1:5000 IF 1:50-200

Anti-TIEG-1/2 KLF10 Antibody (A03419-1) Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using KLF10/11 Antibody. The picture on the right is blocked with the synthesized peptide.

1 Publications Citing This Product

1. PubMed ID: 32156482, Zheng L,Lu H,Li H,Xu X,Wang D.KLF10 is upregulated in osteoarthritis and inhibits chondrocyte proliferation and migration by upregulating Acvr1 and suppressing inhbb expression.Acta Histochem.2020 Apr;122(3):151528.doi:10.1016/j.acthis.2020.151528.Epub 20

Visit bosterbio.com/anti-tieg-1-2-klf10-antibody-a03419-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TIEG-1/2 KLF10 Antibody

For Research Use Only. Not for use in diagnostic procedures.