

Anti-DARPP-32 (phospho-T75) PPP1R1B Antibody

Catalog Number: A03424T75


Overview

Product Name	Anti-DARPP-32 (phospho-T75) PPP1R1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DARPP-32 (phospho-T75) PPP1R1B Antibody catalog # A03424T75. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UD71


Technical Details

Immunogen	Synthetic phosphopeptide derived from human DARPP-32 around the phosphorylation site of Threonine 75.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-DARPP-32 (phospho-T75) PPP1R1B Antibody (A03424T75) Images



Western blot (WB) analysis of p-DARPP-32 (T75) pAb at 1:1000 dilution Lane1:AML-12 whole cell lysate(40ug) Lane2:C6 whole cell lysate(40ug) Lane3:HepG2 whole cell lysate(40ug) Lane4:SGC7901 whole cell lysate(40ug) Lane5:PC3 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of p-DARPP-32 (T75) pAb in paraffin-embedded human colorectal carcinoma tissue at 1:50.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DARPP-32 (phospho-T75) PPP1R1B Antibody

For Research Use Only. Not for use in diagnostic procedures.