

## Anti-SPT5/SUPT5H Antibody Picoband® Biotin Conjugated

Catalog Number: A03426-2-Biotin

### About SUPT5H

Transcription elongation factor SPT5 is a protein that in humans is encoded by the SUPT5H gene. Enables enzyme binding activity and protein heterodimerization activity. Involved in positive regulation of macroautophagy; regulation of RNA metabolic process; and transcription elongation from RNA polymerase II promoter. Located in nucleoplasm. Part of DSIF complex.

### Overview

Product Name	Anti-SPT5/SUPT5H Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O00267

### Technical Details

Immunogen	E.coli-derived human SPT5/SUPT5H recombinant protein (Position: Q533-A1087).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPT5/SUPT5H Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.