

Anti-SPT5/SUPT5H Antibody Picoband® FITC Conjugated

Catalog Number: A03426-2-FITC

About SUPT5H

Transcription elongation factor SPT5 is a protein that in humans is encoded by the SUPT5H gene. Enables enzyme binding activity and protein heterodimerization activity. Involved in positive regulation of macroautophagy; regulation of RNA metabolic process; and transcription elongation from RNA polymerase II promoter. Located in nucleoplasm. Part of DSIF complex.

Overview

| | |
|----------------------|--|
| Product Name | Anti-SPT5/SUPT5H Antibody Picoband® FITC Conjugated |
| Reactive Species | Human, Mouse |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | O00267 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human SPT5/SUPT5H recombinant protein (Position: Q533-A1087). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPT5/SUPT5H Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.