

Anti-SPT5/SUPT5H Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A03426-2-Fluoro488

About SUPT5H

Transcription elongation factor SPT5 is a protein that in humans is encoded by the SUPT5H gene. Enables enzyme binding activity and protein heterodimerization activity. Involved in positive regulation of macroautophagy; regulation of RNA metabolic process; and transcription elongation from RNA polymerase II promoter. Located in nucleoplasm. Part of DSIF complex.

Overview

Product Name	Anti-SPT5/SUPT5H Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O00267

Technical Details

Immunogen	E.coli-derived human SPT5/SUPT5H recombinant protein (Position: Q533-A1087).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPT5/SUPT5H Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.