

Anti-TROVE2/SS-A/RO60 Antibody Picoband® Biotin Conjugated

Catalog Number: A03428-1-Biotin

About RO60

Enables RNA binding activity. Involved in cellular response to interferon-alpha and regulation of gene expression. Located in cytosol and nucleoplasm. Part of ribonucleoprotein complex. The cDNA that encodes the 60-kD Ro/SSA autoantigen recognized by autoantibodies in certain lupus erythematosus (152700) and Sjogren syndrome (270150) sera was cloned and sequenced by Deutscher et al. (1988) and confirmed with minor differences by Ben-Chetrit et al. (1989).

Overview

Product Name	Anti-TROVE2/SS-A/RO60 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P10155

Technical Details

Immunogen	E.coli-derived human TROVE2/SS-A/RO60 recombinant protein (Position: K88-M419).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TROVE2/SS-A/RO60 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.