

## Anti-Eotaxin-3 CCL26 Antibody

Catalog Number: A03434-1

### Overview

Product Name	Anti-Eotaxin-3 CCL26 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Eotaxin-3 CCL26 Antibody catalog # A03434-1. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y258

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CCL26. AA range:31-80
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC-P 1:100-300 ELISA 1:20000

## Anti-Eotaxin-3 CCL26 Antibody (A03434-1) Images




Figure 2. Immunohistochemistry validation of CCL26 using Anti-Eotaxin-3 CCL26 Antibody (A03434-1).

Immunohistochemistry (IHC) analysis of paraffin-embedded Rat Brain




Figure 3. Immunohistochemistry validation of CCL26 using Anti-Eotaxin-3 CCL26 Antibody (A03434-1).

Immunohistochemistry (IHC) analysis of paraffin-embedded Rat Brain




Figure 4. Immunohistochemistry validation of CCL26 using Anti-Eotaxin-3 CCL26 Antibody (A03434-1).

Immunohistochemistry (IHC) analysis of paraffin-embedded Rat Brain




Figure 5. Immunohistochemistry validation of CCL26 using Anti-Eotaxin-3 CCL26 Antibody (A03434-1).

Immunohistochemistry (IHC) analysis of paraffin-embedded Mouse Brain




Figure 6. Immunohistochemistry validation of CCL26 using Anti-Eotaxin-3 CCL26 Antibody (A03434-1).

Immunohistochemistry (IHC) analysis of paraffin-embedded Mouse Brain




Figure 1. Western blotting validation for Anti-Eotaxin-3 CCL26 Antibody A03434-1

Western Blot (WB) analysis of PC-3

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Eotaxin-3 CCL26 Antibody