

Anti-RAGE/AGER Antibody Picoband® Biotin Conjugated

Catalog Number: A03438-1-Biotin

About AGER

The receptor for advanced glycation end products (RAGE) is a multi-ligand member of the immunoglobulin superfamily of cell surface molecules. It interacts with distinct molecules implicated in homeostasis, development and inflammation, and certain diseases such as diabetes and Alzheimer's disease. RAGE is also a central cell surface receptor for amphoterin and EN-RAGE. And RAGE is associated with sustained NF-kappaB activation in the diabetic microenvironment and has a central role in sensory neuronal dysfunction. Moreover, RAGE propagates cellular dysfunction in several inflammatory disorders and diabetes, and it also functions as an endothelial adhesion receptor promoting leukocyte recruitment.

Overview

Product Name	Anti-RAGE/AGER Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q15109

Technical Details

Immunogen	E.coli-derived human RAGE/AGER recombinant protein (Position: K43-P404).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

2 Publications Citing This Product

1. PubMed ID: 10.1016/j.jep.2013.11.015, The anti-inflammation effect of Moutan Cortex on advanced glycation end products-induced rat mesangial cells dysfunction and High-glucose-fat diet and streptozotocin-induced diabetic nephropathy rats
2. PubMed ID: 23219582, Feng L, Zhu M, Zhang M, Jia X, Cheng X, Ding S, Zhu Q. Int Immunopharmacol. 2013 Feb;15(2):206-16. Doi: 10.1016/J.Intimp.2012.11.015. Epub 2012 Dec 5. Amelioration Of Compound 4,4'-Diphenylmethane-Bis(Methyl)Carbamate On High Mobility Group Box1-M...

Visit bosterbio.com/anti-rage-ager-picoband-trade-antibody-a03438-1-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAGE/AGER Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.