

## Anti-Prosurfactant Protein B/SFTPB Antibody Picoband® Biotin Conjugated

Catalog Number: A03441-1-Biotin

### About SFTPB

Pulmonary surfactant-associated protein B is a protein that in humans is encoded by the SFTPB gene. This gene encodes the pulmonary-associated surfactant protein B (SPB), an amphipathic surfactant protein essential for lung function and homeostasis after birth. Pulmonary surfactant is a surface-active lipoprotein complex composed of 90% lipids and 10% proteins which include plasma proteins and apolipoproteins SPA, SPB, SPC and SPD. The surfactant is secreted by the alveolar cells of the lung and maintains the stability of pulmonary tissue by reducing the surface tension of fluids that coat the lung. The SPB enhances the rate of spreading and increases the stability of surfactant monolayers in vitro. Multiple mutations in this gene have been identified, which cause pulmonary surfactant metabolism dysfunction type 1, also called pulmonary alveolar proteinosis due to surfactant protein B deficiency, and are associated with fatal respiratory distress in the neonatal period. Alternatively spliced transcript variants encoding the same protein have been identified.

### Overview

Product Name	Anti-Prosurfactant Protein B/SFTPB Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P07988

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SFTPB, which shares 80% amino acid (aa) sequence identity with both mouse and rat SFTPB.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.

Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: PMID:26097528, 3,4-dihydroxyphenylethanol suppresses irradiation-induced pulmonary fibrosis in adult rats

Visit [bosterbio.com/anti-sftpb-picoband-trade-antibody-a03441-1-boster.html](https://bosterbio.com/anti-sftpb-picoband-trade-antibody-a03441-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prosurfactant Protein B/SFTPB Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.