

## Anti-HSPA2 Antibody

Catalog Number: A03474


### Overview

Product Name	Anti-HSPA2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSPA2 Antibody catalog # A03474. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P54652


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 IHC: 1:50-1:200

## Anti-HSPA2 Antibody (A03474) Images



Western blot analysis of HSPA2 on different lysates using anti-HSPA2 antibody at 1/1,000 dilution. Positive control:  
Lane 1: MCF-7 Lane 2: Hela Lane 3: A431



Immunohistochemical analysis of paraffin-embedded mouse colon tissue using anti-HSPA2 antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-HSPA2 Antibody

For Research Use Only. Not for use in diagnostic procedures.