

Anti-Casein Kinase II Beta (phospho-S209) CSNK2B Antibody

Catalog Number: A03475S209


Overview

Product Name	Anti-Casein Kinase II Beta (phospho-S209) CSNK2B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Casein Kinase II Beta (phospho-S209) CSNK2B Antibody catalog # A03475S209. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P67870

Technical Details

Immunogen	Synthetic phosphopeptide derived from human Cytokeratin 8 around the phosphorylation site of Serine 432
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 96% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Casein Kinase II Beta (phospho-S209) CSNK2B Antibody (A03475S209) Images



Western blot (WB) analysis of Casein Kinase Iibeta (phospho-S209) polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate(40µg) Lane2:Panc1 whole cell lysate(40µg) Lane3:U-87MG whole cell lysate(40µg) Lane4:The Brain tissue lysate of Rat(40µg) Lane5:The Kidney tissue lysate of Mouse(40µg)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Casein Kinase II Beta (phospho-S209) CSNK2B Antibody

For Research Use Only. Not for use in diagnostic procedures.