

## Anti-Desmocollin-1 DSC1 Antibody

Catalog Number: A03478

### About DSC1

Myosins are actin-based motor molecules with ATPase activity. Unconventional myosins serve in intracellular movements. May be involved in the remodeling of the actin cytoskeleton. Binds actin with high affinity both in the absence and presence of ATP and its mechanochemical activity is inhibited by calcium ions. Also acts as a GTPase activating protein on Rho.

Wirth J.A., J. Cell Sci. 109:653-661(1996).

Grewal P.K., Gene 240:389-398(1999).

Post P.L., J. Cell Sci. 111:941-950(1998).

### Overview

Product Name	Anti-Desmocollin-1 DSC1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Desmocollin-1 DSC1 Antibody catalog # A03478. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q08554

### Technical Details

Immunogen	Synthesized peptide derived from human Casein Kinase IIbeta around the phosphorylation site of S209.
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-

	specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows:</p> <p>Boster Bio's internal QC testing used:</p> <p>WB: 1:500-1:2000</p>

## Anti-Desmocollin-1 DSC1 Antibody (A03478) Images




Figure 1. Western blotting validation for Anti-Desmocollin-1 DSC1 Antibody A03478

Western blot analysis of extracts of various cells

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Desmocollin-1 DSC1 Antibody