

## Anti-MKP-2 DUSP4 Antibody

Catalog Number: A03486


### Overview

Product Name	Anti-MKP-2 DUSP4 Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-MKP-2 DUSP4 Antibody catalog # A03486. Tested in WB applications. This antibody reacts with Human, Monkey, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13115

### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human DUSP4. AA range:81-130
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	WB, 1:500-1:2000

## Anti-MKP-2 DUSP4 Antibody (A03486) Images



Western Blot (WB) analysis of HepG2-UV cells using MKP-2 Polyclonal antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MKP-2 DUSP4 Antibody