

Anti-Sulfotransferase 1A1 SULT1A1 Antibody

Catalog Number: A03500


Overview

Product Name	Anti-Sulfotransferase 1A1 SULT1A1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Sulfotransferase 1A1 SULT1A1 Antibody catalog # A03500. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50225

Technical Details

Immunogen	Recombinant fusion protein of human SULT1A1(NP_803565.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Sulfotransferase 1A1 SULT1A1 Antibody (A03500) Images



Western blot analysis of SULT1A1 expression in mouse liver (A), rat liver (B) whole cell lysates.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Sulfotransferase 1A1 SULT1A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.