

Anti-Eotaxin-2 CCL24 Antibody

Catalog Number: A03513


Overview

Product Name	Anti-Eotaxin-2 CCL24 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Eotaxin-2 CCL24 Antibody catalog # A03513. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00175


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 61-110
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

Anti-Eotaxin-2 CCL24 Antibody (A03513) Images



Immunohistochemistry validation of CCL24 using Anti-Eotaxin-2 CCL24 Antibody (A03513). Immunohistochemical analysis of paraffin-embedded human-lung



Immunohistochemistry validation of CCL24 using Anti-Eotaxin-2 CCL24 Antibody (A03513). Immunohistochemical analysis of paraffin-embedded human-lung

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Eotaxin-2 CCL24 Antibody

For Research Use Only. Not for use in diagnostic procedures.