

## Anti-Synuclein-gamma SNCG Monoclonal Antibody

Catalog Number: A03523


### Overview

Product Name	Anti-Synuclein-gamma SNCG Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Synuclein-gamma SNCG Monoclonal Antibody catalog # A03523. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 1H1D2
Formulation	Purified antibody in PBS containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	O76070


### Technical Details

Immunogen	Purified recombinant fragment of Synuclein-gamma expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

## Anti-Synuclein-gamma SNCG Monoclonal Antibody (A03523) Images



Western Blot analysis using Synuclein-gamma Monoclonal Antibody against truncated Synuclein-gamma recombinant protein.



Immunohistochemistry analysis of paraffin-embedded human ovary carcinoma (left) and breast carcinoma (right), showing cytoplasmic (ovary carcinoma) localization, cytoplasmic and nuclear (breast carcinoma) localization with DAB staining using Synuclein-gamma Mo

### 1 Publications Citing This Product

1. PubMed ID: -, Ou Qiao,Xinyu Zhang,Yi Zhang et al.Cerebralcare Granule® enhances memantine hydrochloride efficacy in APP/PS1 mice by ameliorating amyloid pathology and cognitive functions,08 April 2021,PREPRINT (Version 1) available at Research Square [https://doi.org/10.21203/rs.3.rs-366097/v1]

Visit [bosterbio.com/anti-synuclein-gamma-antibody-a03523-boster.html](https://bosterbio.com/anti-synuclein-gamma-antibody-a03523-boster.html) to see all 1 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Synuclein-gamma SNCG Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.