

Anti-Synoretin (A111) SNCG Antibody

Catalog Number: A03523-1


Overview

Product Name	Anti-Synoretin (A111) SNCG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Synoretin (A111) SNCG Antibody catalog # A03523-1. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O76070

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 81-130 of Human Synoretin.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-Synoretin (A111) SNCG Antibody (A03523-1) Images



Immunohistochemistry (IHC) analyzes of Synoretin (A111) pAb in paraffin-embedded human lung carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Synoretin (A111) SNCG Antibody

For Research Use Only. Not for use in diagnostic procedures.