

Anti-Synoretin (A111) SNCG Antibody

Catalog Number: A03523-1

About SNCG

Displays serine/threonine-specific phosphorylation of myelin basic protein and histone (MBP) in vitro.

Okano I., J. Biol. Chem. 270:31321-31330(1995).

Osada H., Biochem. Biophys. Res. Commun. 229:582-589(1996).

Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).

Overview

Product Name	Anti-Synoretin (A111) SNCG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Synoretin (A111) SNCG Antibody catalog # A03523-1. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O76070

Technical Details

Immunogen	Synthesized peptide derived from human PSYR protein.
Predicted Reactive Species	Boar, Bovine, Canine, Golden Hamster
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
IHC: 1:50-1:200

Anti-Synoretin (A111) SNCG Antibody (A03523-1) Images




Figure 1. Immunohistochemistry validation of SNCG using Anti-Synoretin (A111) SNCG Antibody (A03523-1).

Immunohistochemistry (IHC) analyzes of Synoretin (A111) pAb in paraffin-embedded human lung carcinoma tissue at 1:100.

For more protocol information of IHC

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Synoretin (A111) SNCG Antibody