

Anti-I10R2 IL10RB Antibody

Catalog Number: A03536

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-I10R2 IL10RB Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-I10R2 IL10RB Antibody catalog # A03536. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q08334

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 240-320
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	I10R2 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000 RRM2B/p53R2 antibody is tested in Western Blot and ELISA and suitable for IP, IHC, and IF. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 41 kDa in size corresponding to RRM2B1 by western blotting in the appropriate cell lysate or extract.

Anti-IL10R2 IL10RB Antibody (A03536) Images



Anti-IL10R2 IL10RB Antibody

For Research Use Only. Not for use in diagnostic procedures.