

Anti-40S ribosomal protein S2 RPS2 Antibody

Catalog Number: A03548-2


Overview

| | |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-40S ribosomal protein S2 RPS2 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-40S ribosomal protein S2 RPS2 Antibody catalog # A03548-2. Tested in WB,ICC,IP applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IP, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P15880 |

Technical Details

| | |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen | Synthetic peptide, corresponding to amino acids 248-293 of Human RPS2. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 ICC: 1:50-1:200 IP: 1:50-1:200 |

Anti-40S ribosomal protein S2 RPS2 Antibody (A03548-2) Images



Western blot (WB) analysis of RPS2 polyclonal antibody at 1:500 dilution Lane1:HEK293T cell lysate Lane2:sp2/0 cell lysate Lane3:H9C2 cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-40S ribosomal protein S2 RPS2 Antibody

For Research Use Only. Not for use in diagnostic procedures.