

## Anti-NHEJ1 Antibody

Catalog Number: A03552

### Overview

Product Name	Anti-NHEJ1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NHEJ1 Antibody catalog # A03552. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H9Q4

### Technical Details

Immunogen	Recombinant fusion protein of human NHEJ1(NP_079058.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:200-1:2000

## Anti-NHEJ1 Antibody (A03552) Images




Figure 1. Western blotting validation for Anti-NHEJ1 Antibody A03552

Western blot analysis of extracts of various cells

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NHEJ1 Antibody