

## Anti-Glycine N-methyltransferase GNMT Antibody

Catalog Number: A03562-1


### Overview

Product Name	Anti-Glycine N-methyltransferase GNMT Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-Glycine N-methyltransferase GNMT Antibody catalog # A03562-1. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14749

### Technical Details

Immunogen	Recombinant fusion protein of human GNMT(NP_061833.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-Glycine N-methyltransferase GNMT Antibody (A03562-1) Images



Western blot analysis of extracts of various cell lines, using GNMT antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glycine N-methyltransferase GNMT Antibody

For Research Use Only. Not for use in diagnostic procedures.