

## Anti-TAF II p68 (N380) TAF15 Antibody

Catalog Number: A03567


### Overview

Product Name	Anti-TAF II p68 (N380) TAF15 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TAF II p68 (N380) TAF15 Antibody catalog # A03567. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92804


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 347-398 of Human TAF II p68.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-TAF II p68 (N380) TAF15 Antibody (A03567) Images



Western blot (WB) analysis of TAF II p68 (N380) pAb at 1:500 dilution Lane1:3T3-L1 whole cell lysate(40ug) Lane2:C6 whole cell lysate(40ug) Lane3:MCF-7 whole cell lysate(40ug) Lane4:HEK293T whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of TAF II p68 (N380) pAb in paraffin-embedded human tonsil carcinoma tissue at 1:50.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TAF II p68 (N380) TAF15 Antibody

For Research Use Only. Not for use in diagnostic procedures.