

## Anti-Dishevelled 3/DVL3 Antibody Picoband® Biotin Conjugated

Catalog Number: A03577-Biotin

### About DVL3

Segment polarity protein dishevelled homolog DVL-3 is a protein that in humans is encoded by the DVL3 gene. It is mapped to 3q27.1. This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.

### Overview

Product Name	Anti-Dishevelled 3/DVL3 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q92997

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Dishevelled 3, identical to the related mouse sequence.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Dishevelled 3/DVL3 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.