

## Anti-Dishevelled 3/DVL3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A03577-Fluoro488

### About DVL3

Segment polarity protein dishevelled homolog DVL-3 is a protein that in humans is encoded by the DVL3 gene. It is mapped to 3q27.1. This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.

### Overview

Product Name	Anti-Dishevelled 3/DVL3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q92997

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human Dishevelled 3, identical to the related mouse sequence.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Dishevelled 3/DVL3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.