

Anti-GRF-1 (phospho-Y1087) ARHGAP35 Antibody

Catalog Number: A03592Y1087


Overview

Product Name	Anti-GRF-1 (phospho-Y1087) ARHGAP35 Antibody
Reactive Species	Human
Description	Boster Bio Anti-GRF-1 (phospho-Y1087) ARHGAP35 Antibody catalog # A03592Y1087. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NRY4

Technical Details

Immunogen	Synthetic phosphopeptide derived from human GRF-1 around the phosphorylation site of Tyrosine 1087.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-GRF-1 (phospho-Y1087) ARHGAP35 Antibody (A03592Y1087) Images



Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GRF-1 (phospho-Y1087) ARHGAP35 Antibody

For Research Use Only. Not for use in diagnostic procedures.