

Anti-Cytochrome P450 26A1 CYP26A1 Antibody

Catalog Number: A03646


Overview

Product Name	Anti-Cytochrome P450 26A1 CYP26A1 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Cytochrome P450 26A1 CYP26A1 Antibody catalog # A03646. Tested in WB,IHC applications. This antibody reacts with Human,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43174


Technical Details

Immunogen	Recombinant protein corresponding to human CYP26A1 aa 200-400.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200

Anti-Cytochrome P450 26A1 CYP26A1 Antibody (A03646) Images



Western blot analysis of CYP26A1 on HepG2 cell lysate using anti-CYP26A1 antibody at 1/500 dilution.



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-CYP26A1 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 26A1 CYP26A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.