

Anti-Cytochrome P450 26A1 CYP26A1 Antibody

Catalog Number: A03646-1

Overview

| | |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Anti-Cytochrome P450 26A1 CYP26A1 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Cytochrome P450 26A1 CYP26A1 Antibody catalog # A03646-1. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat. |
| Application | ELISA, IF, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | O43174 |

Technical Details

| | |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Immunogen | The antiserum was produced against synthesized peptide derived from human Cytochrome P450 26A1. AA range:251-300 |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:20000 |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 26A1 CYP26A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.