

Anti-ULBP2 Antibody

Catalog Number: A03657-1

About ULBP2

NKG2D is an activating cell surface receptor that is predominantly expressed on cytotoxic immune cells. NKG2D recognizes a wide range of ligands. In humans, the NKG2D ligand (NKG2DL) includes MICA, MICB, and six members of the ULBP family (ULBP1-6) (PMID: 29568297; 30813924). ULBPs are MHC class I-related molecules (PMID: 11239445). ULBP2 contains two Ig-like domains, the alpha-1 and alpha-2 domains characteristic of the MHC class I family, but lacks the alpha-3 domain. It can be expressed as either a transmembrane or GPI-linked protein, or released from the cell surface (PMID: 21224393; 24614922). ULBP2 binds and activates the NKG2D, mediating the recruitment and activation of NK cells.

Overview

Product Name	Anti-ULBP2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-ULBP2 Antibody catalog # A03657-1. Tested in WB, IHC, ELISA applications. This antibody reacts with Human.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BZM5

Technical Details

Immunogen	E.coli-derived human ULBP2 recombinant protein (Position: R27-Q221).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ULBP2 Antibody

For Research Use Only. Not for use in diagnostic procedures.